



Site Preparation Specification

Purpose of Procedure

To ensure that the installation site is properly evaluated and prepared with the appropriate utilities, consumables and supplies for the successful installation of AGILENT instruments and systems.

Procedure Checklist



AGILENT G2445 Mainframe:

footprint:

Weight:	80 kg	Height:	55.2 cm
	176 lb		21.7 in
Depth:	69.6 cm	Width:	75.4 cm
	27.4 in		29.7 in

Tick Boxes

Customer Responsibilities

Customers should ensure that all necessary operating supplies, consumables and usage dependent items such as columns, vials, syringes and solvents required for the successful installation of instruments and systems are available. Installation sites should be prepared in accordance with the following specifications. An AGILENT customer engineer will call approximately 2 weeks prior to installation to confirm site readiness.

Maximum cabinet dimensions:*

Weight:	83.0 kg	Height:	64.2 cm
	182.0 lb		25.2 in
Depth:	75.1 cm	Width:	82.0 cm
	29.6 in		32.2 in

E1M18 Mechanical Pump:

Weight:	32.0 kg	Height:	23 cm
	70.4 lb		9.2 in
Depth:	51.0 cm	Width:	17.0 cm
	20.4 in		6.8 in

AGILENT G1947A APCI Interface:**

Weight:	1.7 kg	Height:	23 cm
	3.75 lb		9.2 in
Depth:	13.0 cm	Width:	18 cm
	5.1 in		7.1 in

Important Information

This checklist is designed to be used in conjunction with the AGILENT 1100 Series LC/MSD Site Preparation Manual, G2440-90004. If you have problems providing any of the following, please contact your local AGILENT office for assistance. Assistance with user specific applications may be provided but should be contracted separately. Users of the instrument should be present throughout the installation and familiarization otherwise important operational, maintenance and safety information may be missed.

AGILENT G1948A API-ES Interface:**


Weight:	1.7 kg	Height:	17 cm
	3.75 lb		6.8 in
Depth:	9.5 cm	Width:	18.0 cm
	3.7 in		7.1 in

* Maximum cabinet dimensions are for an LC/MSD Trap G2440AA with G1947A APCI interface installed.

** At least 30 cm (1 ft) to the left of the instrument and at least 5 cm (2 in) to the rear should be added to the maximum cabinet dimensions to provide adequate airflow and instrument access.

Site Preparation Specification

Environmental Conditions



Temperature: 15 to 35 °C (59 to 95 °F)
at constant temperature (variations <3 °C/hr).
Analytical specifications will be met within the temperature
range of 21 °C ± 3 °C (70 °F ± 6 °F)
Humidity: <95% relative, non-condensing.



Power



U.S. & Japan: 200 - 210 VAC, +5%/-10%
50/60 Hz ± 5%
2000 VA max
Europe: 220 - 240 VAC, +5%/-10%
50/60 Hz ± 5%
2000 VA max




Heat Dissipation



2000 Watts (6800 BTU / hour)



Nitrogen Gas Supply



Purity: 99.5% or better - Cylinder
98.0% or better - Nitrogen gas generator or liquid nitrogen
Dewar
The balance should consist of oxygen and/or argon and must
be hydrocarbon free (<0.1 ppm hydrocarbons).
Pressure: 80-100 psi. One 1/4" Swagelok.
Volume: Up to 15 liters/min. outlet fitting is required to
connect the LC/MSD



Helium Gas Supply

Purity: 99.99% or better
Pressure: 60-80 psi maximum
Helium Usage: < 5sccm/min (standard cubic
Centimeter)



Exhaust Venting Requirements

Capacity: Up to 15 liters/min. with separate exhaust lines for
pump and API source.



Remote Diagnostics

- One analog phone line is recommended to provide
remote diagnostics capability for the LC/MSD Trap



Site Preparation Specification

- A second phone line is also strongly recommended for communication with the system operator



Solvents and Supplies

- Iso-Propanol, Methanol and Acetonitrile: HPLC grade purity is required
- MiliQ water or HPLC Grade water is recommended
- Buffer: Acetic acid (recommended), Ammonium

Format

- Eppendorf pipettes (10-100µl and 100-1000µl)
- Eppendorf pipette tip

Note:

This table lists operation power for G2440AA LC/MSD Trap and its accessories:

Line Voltage	120V AC	230V AC
Mass Spectrometer	200V - 210V AC +5% - 10% 50/60Hz ± 5%	220V - 240V AC +5% - 10% 50/60Hz ± 5%
	Dual Phase 16A Max	Single Phase 16A Max
Computer Hardware	100V - 240V AC 50/60 Hz (110V AC)	100V - 240V AC 50/60 Hz (230V AC)
	5A, 2 outlets	3A, 2 outlets
Printer LaserJet (Purchase separately)	100V - 120V AC 50/60 Hz	220V - 240V AC 50 Hz
	3.3A, 1 outlet	1.7A, 1 outlet
Syringe Pump (Include with G2440AA)	115VAC – 230VAC 50/60 Hz 0.1A, 1 outlet	
LC 1100 (Purchase separately)	100V - 240V AC 50/60 Hz	
Number of outlets varies with actual configuration		